
Vicks® Vaporizer

Owner's Manual

PARTS DIAGRAM:

- | | |
|---------------------|-----------------------------|
| A. Finger Grips | E. In Canada (Medicine Cup) |
| B. Steam Guard | In U.S. (Water Container) |
| C. Outer Cap | F. Protective Housing |
| D. Night Light Lens | |

IMPORTANT! For the most effective and safe use of our vaporizer, please **READ AND SAVE THESE INSTRUCTIONS...** especially those regarding the addition of salt to the water. In "soft water" areas, salt must be added to water to generate steam. In "hard water" areas, the vaporizer may operate without the addition of salt.

WARNING: KEEP OUT OF REACH OF CHILDREN. STEAM CAN CAUSE SERIOUS BURNS.

We recognize that many of our customers use vaporizers in homes with young children. Carefully supervise your children when using a vaporizer, especially crawling infants and toddlers. Please be sure to take the time to instruct children that a vaporizer is not a toy. Vaporizers produce hot steam and could cause severe burns and injuries.

Vaporizer should always be placed on a firm, flat, waterproof surface at least four feet away from bedside and out of reach of patient and children. Be sure the vaporizer is in a stable position and power cord is out of the way to prevent vaporizer from being upset or tipped over. Protect floor, rug or furniture by placing the vaporizer on a non-metallic, waterproof, heat-resistant surface.

ADDING SALT FOR STEAM FLOW: Steam output depends directly on the mineral content of the water and the amount of salt you add. For the vaporizer to operate in areas of low or moderate mineral-content water, 3 TO 4 PINCHES OF SALT (1/8 teaspoon) MUST BE ADDED TO THE WATER AND MIXED THOROUGHLY.

If an adequate steam flow is not obtained after ten minutes, add another PINCH or two as needed, mix thoroughly. **CAUTION:** Add salt only by the PINCH, since too much salt can cause excessive boiling, a blown fuse or tripped circuit breaker.

If you live in a "hard-water" area, the vaporizer may operate well without salt additive. The salt and other minerals in your water remain inside the vaporizer and do not affect the purity of the steam.

USE FINGER GRIPS when handling steam unit, grasping from the rear of the outer cap. See Figure C. Do not move vaporizer while in use.

If you have any questions about the operation of your vaporizer, call our toll-free Consumer Relations Line: 1-800-477-0457.

INSTRUCTIONS:

- UNRAVEL CORD BEFORE USE.** Failure to do so could cause overheating and possible fire hazard.

Always unplug cord from electrical outlet before removing steam unit or when not in use. Allow vaporizer to cool before handling.

- FILL WATER CONTAINER.** Remove TWIST AND LIFT steam unit by aligning arrows on the head and the bowl. (Figure A) Fill container with tap water to WATER LEVEL MARKING on bowl. (Figure B) DO NOT OVERFILL. Vicks Vaporizer may also be filled in a large sink or bathtub.

- ADD SALT.** For tap water of low or moderate mineral content, ADD 4 TO 5 PINCHES (1/8 teaspoon) of ordinary table salt to the water and mix thoroughly. Salt has no effect on the purity of the steam or medicine. **CAUTION:** Too much salt can cause excessive boiling or blown fuse.

- REPLACE STEAM UNIT.** Hold down steam unit for ten seconds to allow water to enter the housing and seat it properly in water container. If steam unit floats, container is overfilled. (Figure C)

- PLACE VAPORIZER IN A SAFE POSITION.** Vaporizer should always be placed on a firm, flat surface at least four feet away from bedside and out of the reach of person and children. Be sure the vaporizer is in a stable position and power cord is out of the way to prevent vaporizer from being upset or tipped over.

Protect floor, rug or furniture by placing the vaporizer on a non-metallic, waterproof, heat-resistant surface.

CAUTION: Do not direct steam flow directly at person, walls or furniture. You will still obtain desired level of humidity. **KEEP OUT OF REACH OF CHILDREN. STEAM CAN CAUSE SERIOUS BURNS.**

- VAPORIZER INHALANTS:**

In U.S. – For best results, fill the medicine cup on the steam unit (located directly above the Steam Guard grill) with one tablespoon of Kaz Inhalant. Vicks VapoSteam® should be put directly into the vaporizer water container (see VapoSteam package directions). (Figure D)

In Canada – For best results, fill the medicine cup on the steam unit (located directly above the Steam Guard grill) with one tablespoon of Vicks VapoSteam® (see VapoSteam package directions). (Figure D)

If you use Vicks VapoSteam®, follow directions on the VapoSteam package. Please pay careful attention to the added cleaning directions in the *Daily Maintenance and Care and Cleaning* section.

- PLUG IN.** Plug power cord into 110-120 Volts AC outlet in good condition, fused for at least 15 Amps. Do not use with extension cords.

This product has a polarized plug (one blade is wider than the other) as a safety feature. This plug will fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. **DO NOT ATTEMPT TO DEFEAT THIS SAFETY FEATURE.**

Vaporizer will produce a gentle flow of steam within a few minutes. If steaming does not start within ten minutes, you may need to add salt. Please see the Troubleshooting checklist. The vaporizer will automatically stop producing when water has reached minimum level (about 3/4 inch from bottom of container). The night-light will remain on, indicating that the vaporizer is plugged in. ALWAYS unplug steam unit and allow it to cool before removing from water container.

Do not lift or move vaporizer while it is in use.

- USE FRESH WATER AND FRESH INHALANT** each time vaporizer is used.
- AFTER EACH USE.** Empty any unused water from the water container. Shake container gently back and forth to remove any remaining water. Rinse container thoroughly, drain and wipe with a dry, clean cloth.

NOTE: If you are using Vicks VapoSteam be sure any residue is thoroughly cleaned. Refer to *Care and Cleaning* section.

IMPORTANT: To maintain proper operation of your vaporizer, be sure to follow *Care and Cleaning* instructions.

SPECIAL NOTES

Discoloration/Black Particles

After use, the water residue in the bottom of the container may contain black particles, and, in some areas the water and the bottom of container may discolor. These conditions are normal and do not effect the purity of the steam or medicine. The container is easily cleaned with soap and water. Rinse and dry the container after cleaning.

Ventilation

When vaporizer is used in closed or small rooms, check periodically for proper ventilation, as excessive humidity can cause condensation stains on walls.

Warning

KEEP OUT OF REACH OF CHILDREN. Carefully supervise children when the vaporizer is operating, especially infants and toddlers. Vaporizers produce hot steam that can cause serious burns.

Condition of Sale

As a condition of sale, the purchaser assumes responsibility for the proper care and use of this appliance in accordance with our printed instructions. The purchaser or user must judge for himself or herself when to use it and length of use.

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
TROUBLESHOOTING: WHAT TO DO IF...

A. Vaporizer does not operate to steams too slowly:

1. Make sure vaporizer is plugged in and the night-light is lit. If the night-light is not on, you may have a blown fuse, tripped circuit breaker or a defective electrical outlet. Try a different outlet. You may want to contact an electrician to check the potentially defective outlet.
2. The fault may be the low mineral content of the water in your area. Unplug vaporizer from electrical outlet and allow it to cool. Twist and Lift to remove steam unit slowly to allow hot water to first drain out. Add one or two PINCHES of salt to the water and mix thoroughly. Replace the steam unit and plug in vaporizer. It should steam within 5 minutes. Add salt only by the PINCH, since too much salt can cause excessive boiling or blowing of fuse.
3. If vaporizer has been operating normally and is now failing to steam, it needs to be cleaned. See the *Weekly Maintenance* section of *Care and Cleaning* instructions.

B. Vaporizer operates too fast, fuse blows:

Overboiling, flickering lights, blown fuse or heavy steam flow generally indicates one of the following causes...

1. Too much salt was added to the water.

To correct the condition; Disconnect vaporizer from electrical outlet and allow it to cool. Pour out water and rinse container. Rinse steam unit as described in the *Weekly Maintenance* section of *Care and Cleaning* instructions. Refill with fresh water but **DO NOT ADD SALT**. If vaporizer now does not operate or steams too slowly, see above directions in section A.

2. You live in a "hard water" area and minerals in the tap water are causing the vaporizers to steam too rapidly.

To correct the condition: Disconnect vaporizer from electrical outlet and allow it to cool. Remove steam unit allowing hot water to first dribble out, pour out water and rinse container. Rinse steam unit as described in the *Weekly Maintenance* section of *Care and Cleaning* instructions. Refill with fresh water but **DO NOT ADD SALT**.

3. If unit still steams too rapidly, use one-half distilled water and one-half tap water. Food and drug stores usually stock distilled water.

In "hard water" areas it is advisable to clean the unit every few days in accordance with the *Weekly Maintenance* section of *Care and Cleaning* instructions.

C. Vaporizer is making crackling/buzzing/humming sounds or water container is unusually hot:

1. You may have added too much salt. See section B above.
2. You may live in a "hard water" area. See section B above.
3. This condition can happen when the water level in the container is low. Disconnect vaporizer from electrical outlet, allow to cool. Twist and Lift to remove steam unit slowly to allow hot water to first drain out. Fill container with tap water to water level marking.
4. This may be an indication your vaporizer needs cleaning. Refer to the *Weekly Maintenance* section of *Care and Cleaning* instructions. If you are still experiencing difficulty in operating your vaporizer, please call our Consumer Service Department at 800-477-0457.

LIMITED WARRANTY

This vaporizer is warranted against defective materials and workmanship. See carton for warranty period. Any defective part will be replaced at no charge if it has not been tampered with and appliance has been used according to our printed directions. Warranty gives you specific legal rights: you may have other rights which vary from state to state. If you experience unsatisfactory operations, first follow *Troubleshooting* and *Care and Cleaning* instructions. If the unit still does not operate satisfactorily, call Consumer Relations Department at 1-800-477-0457 or write to Kaz, Inc., Consumer Relations Department, 250 Turnpike Rd., Southborough, MA 01772, explaining just how the vaporizer is operating. Specify model number. Our consumer service representative will provide further instructions on how to correct the trouble yourself, or will ask you to return the steam unit for repair or replacement.

CARE AND CLEANING

ALWAYS DISCONNECT THE PLUG FROM THE ELECTRICAL OUTLET AND ALLOW UNIT TO COOL BEFORE CLEANING. REMOVE STEAM UNIT, LETTING ALL WATER DRAIN BACK INTO THE CONTAINER.

DO NOT ATTEMPT TO DISASSEMBLE THE ELECTRODE HEATING UNIT. DOING SO WILL VOID OUR WARRANTY AND COULD DAMAGE THE VAPORIZER OR CAUSE PERSONAL INJURY.

Daily Maintenance:

Empty unused water from water container. Rinse container thoroughly with water, drain and wipe with a dry, clean cloth.

Weekly Maintenance:

1. Empty unused water from water container. Rinse container thoroughly with water, drain and wipe with a dry, clean cloth.
2. With steam unit safely unplugged, soak steam unit in 4 inches of white vinegar for 10 minutes (Figure E).

Then, under a sink, with fingers covering the bottom hole of the steam unit, add tap water through the Steam Guard. (Figure F) Shake vigorously from side to side and then let water empty through the bottom hole. (Figure G) Repeat several times until black particles are no longer emitted and vinegar odor is gone. If using Vicks VapoSteam, make sure any residue is cleaned off steam unit using a mild detergent.

3. Rinse area under the Steam Guard to remove any dust or dirt particles. With long toothpick, you can gently poke through the steam outlet and draining hole if they've become obstructed.

4. Disinfecting the Reservoir: Fill reservoir with water and add 1 teaspoon of bleach per gallon of water. Swish to wet all inner surfaces. Let stand for 20 minutes. Then empty reservoir and rinse with tap water until bleach smell is gone. Dry with a clean cloth.

End-Of-Season Maintenance/Storage:

Follow *Weekly Maintenance* procedure at end-of season or when the vaporizer will not be used for several days. Be sure to dry all parts and store in a cool, dry location. **DO NOT STORE UNIT WITH WATER INSIDE.**

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©2007
Vicks is a registered trademark of The Procter & Gamble Company.